

ACE Glass Fiber Filter Discs

ACE Designation	Porosity Maximum Pore Diameter Range (micron)		Corning Equivalent	Kimble Equivalent	Uses
A	145-174	ACE	EC (170-220)	EC (170-220)	Coarse Filtration
B	70-100		—	—	Coarse Filtration
C	25-50		C (40-60)	C (40-60)	Gas Dispersion
D	10-20		M (10-15)	M (10-15)	Extraction
E	4-8		F (4-5.5)	F (4-5.5)	Extraction
VF	2-2.5	Robu	VF (2-2.5)	—	Bacteria Filtration
UF	0.9-1.4		UF (0.9-1.4)	—	Bacteria Filtration