

Matching Ace Stir Bearings to the Appropriate Shafts

Bearing Type	Size	Ace Stir Bearing Codes	Use with Ace Code Stirring Shafts Listed Below
Trubore, Glass	5mm	9524-04, 9527-08	9534-04, 9535-06, 9541-04, 9541-15
Trubore, Glass	6mm	9524-06, 9524-08, 9527-12, 9527-14, 9529	9534-06
Vacuum	9mm	8098, 8099, 8133, 9528	8134, 9530
Trubore, Glass	10mm	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8051, 8053, 8055, 8133, 8066	8068, 8070, 8071, 8073, 9532, 9533
Trubore, Glass	10mm	8041	Complete Assemblies
Trubore, Glass	19mm	8059, 8060, 8061, 8065	8076, 8077, 8078, 8079
Glass, Pressure	10mm	8044	8074, 8075
Trubore, PTFE/Glass	10mm	8066 (Plain), 8066 (Debris Trap)	8074, 8075
Trubore, PTFE/Glass	19mm	8067 (Trubore), 8067 (Debris Trap)	8076
Trubore, PTFE/Glass	28mm	8067 (Trubore), 8067 (Debris Trap)	8080
Mechanical Seal/PTFE	10mm	8050-02	8068, 8073, 8074, 8075, 9532, 9533
Mechanical Seal/PTFE	19mm	8050-06	8076, 8077, 8078
Mechanical Seal/PTFE	28mm	8050-08	8080
Mechanical Seal/PTFE	6mm	13443-06, 13443-08	9534-06, 9534-40
Mechanical Seal/PTFE	8mm	13443-10	9534-06, 9534-40
Mechanical Seal/PTFE	10mm	13443-12	8068, 8070, 8071, 8073, 8074, 8075, 9532, 9533
PTFE	6mm	13445-06, 13445-09	9534-40
PTFE	10mm	13445-30, 13445-32, 13445-34, 13445-36	8075, 8075
PTFE	19mm	13445-46	8076